

pH Buffer Solutions

Available in different pack sizes.

Colored Buffers for laboratory use, pH measurements, traceable to SRM from NIST (@20°C)

Tolerance: $\pm 0.02/20^{\circ}\text{C}$

Product code	Product name	Formulation
CL03.0702	Buffer pH 4.00 (colored)	Citrate-NaCl-NaOH (red)
CL03.0704	Buffer pH 7.00 (colored)	Phosphate (green)
CL03.0710	Buffer pH 10.00 (colored)	H_3BO_3 -KCl-NaOH (blue)

Buffers for laboratory use, pH measurements, traceable to SRM from NIST (@20°C)

Tolerance: $\pm 0.02/20^{\circ}\text{C}$

Product code	Product name	Formulation
CL03.0203	Buffer pH 1.00	KCl-HCl
CL03.0209	Buffer pH 2.00	Citrate-NaCl-HCl
CL03.0210	Buffer pH 3.00	Citrate-NaCl-NaOH
CL03.0213	Buffer pH 4.00	Citrate-NaCl-NaOH
CL03.0214	Buffer pH 5.00	Citrate-NaOH
CL03.0215	Buffer pH 6.00	Citrate-NaOH
CL03.0216	Buffer pH 7.00	Phosphate
CL03.0217	Buffer pH 8.00	H_3BO_3 -KCl-NaOH
CL03.0218	Buffer pH 9.00	H_3BO_3 -KCl-NaOH
CL03.0204	Buffer pH 10.00	H_3BO_3 -KCl-NaOH
CL03.0206	Buffer pH 11.00	H_3BO_3 -KCl-NaOH
CL03.0207	Buffer pH 12.00	Phosphate-NaOH
CL03.0208	Buffer pH 13.00	Glycin-NaCl-NaOH

Contact your
distributor

We have the solutions!